

Labeling experiment with stable isotopes

General considerations for stable isotope (e.g. C-13, N-15, H-2 (D)) labeling:

- The minimum labeling duration depends on the metabolic pathway of interest. For instance, metabolites in the glycolysis pathway only take minutes to reach isotopic steady-state for many cell types, while other metabolic pathways (e.g. lipids) might take days.
- plate about 200k (or more) cells per well in 6well plates (3 wells per condition minimum) and incubate o/n (for adherent cells, add additional well(s) for cell counts)
 - after 24h rinse with 1x PBS and add fresh medium (with heavy tracer or not): 1.5-2 ml is sufficient to cover the cells for 24h, but keep in mind that the cells might deplete nutrients within that amount of time

Preparation of medium (metabolic footprint) and plasma samples

- remember that you will need 3 extra samples with 'fresh' medium: that is unused, spare medium you put on your cells
- take 20 µl medium (or blood serum) and 5 nmol norvaline to 500 µl 80% MeOH, vortex and spin for 5' at top speed
- transfer supernatant into glass tube and dry as described above
- keep samples at -80C at CNSI

Equipment and reagents needed for this protocol

- | | | | |
|--------------------|------------|-------------------|---|
| - Ammonium acetate | A1542-500G | Fisher | for molecular biology, ≥98% |
| - glass vials: | 03-410-151 | Fisher | 1.8 mL Volume; Clear Glass, 12x32 mm, 9 mm thread |
| - caps: | 03-379-123 | Thermo Scientific | Rubber/Silicone Septa |
| - MeOH: | A456-1 | Fisher | Fisher Methanol (Optima* LC/MS) |
| - H2O: | W5-1 | Fisher | Water, Glass Bottle; 1L |
| - Norvaline: | N7502-25G | Sigma | DL-Norvaline |

Alternatively: American Chromatography Supplies

- | | | | |
|----------------|-------------|-----|---|
| - glass vials: | VT009M-1232 | ACS | 1.8 mL Volume; Clear Glass, 12x32 mm, 9 mm thread |
| - caps: | C395E-09SB | ACS | Bonded PTFE/Silicone Septa |
| - caps: | C394-09SB | ACS | Bonded PTFE/Rubber Septa |

C-13- and N-15-labeled metabolites (from [Cambridge Isotope Laboratories](#) if not otherwise stated)

- | | | |
|-----------------|------------|-----|
| - U13C Glucose: | CLM-1396-1 | 1 g |
|-----------------|------------|-----|