

Preparing cells for a metabolomics experiment

General considerations for stable isotope (e.g. C-13, N-15, H-2 (D)) labeling:

- The minimum labeling duration depends on the metabolic pathway of interest. For instance, metabolites in the glycolysis pathway only take minutes to reach isotopic steady-state for many cell types, while other metabolic pathways (e.g. lipids) might take days.
- you can label animals through i.p. or i.v. injection (bolus) or over time through a catheter into the tail vein, for instance
 - make sure to have a positive/input control, i.e. blood plasma, to verify that the tracer has gone into the blood stream

Metabolite extraction from tissue*

- cut off 10-20 mg tissue on ice and homogenize in 1 ml 80% MeOH (-80C) (you might have to rinse the tissue to remove blood if possible)
 - spin down at top speed (@4C) 1-2x for 5 min
 - resuspend cell / tissue pellet in protein lysis buffer and measure protein concentration (see above)
 - transfer fraction of supernatant (usually the equivalent of 1-5 µg protein) into glass tube, add 5 nmol norvaline
- dry samples using EZ-2Elite evaporator at 30C using program 3 (aqueous)
- keep dried samples at -80C at CNSI

Equipment and reagents needed for this protocol

| | | | |
|--------------------|------------|-------------------|---|
| - Ammonium acetate | A1542-500G | Fisher | for molecular biology, ≥98% |
| - glass vials: | 03-410-151 | Fisher | 1.8 mL Volume; Clear Glass, 12x32 mm, 9 mm thread |
| - caps: | 03-379-123 | Thermo Scientific | Rubber/Silicone Septa |
| - MeOH: | A456-1 | Fisher | Fisher Methanol (Optima* LC/MS) |
| - H2O: | W5-1 | Fisher | Water, Glass Bottle; 1L |
| - Norvaline: | N7502-25G | Sigma | DL-Norvaline |

| | | | |
|----------------|----------------------------------|-----|---|
| Alternatively: | American Chromatography Supplies | | |
| - glass vials: | VT009M-1232 | ACS | 1.8 mL Volume; Clear Glass, 12x32 mm, 9 mm thread |
| - caps: | C395E-09SB | ACS | Bonded PTFE/Silicone Septa |
| - caps: | C394-09SB | ACS | Bonded PTFE/Rubber Septa |

C-13- and N-15-labeled metabolites (from [Cambridge Isotope Laboratories](#) if not otherwise stated)

- U13C Glucose: CLM-1396-1 1 g

*Thanks to Marcus Seldin and Margarete Mehrabian